

Post-IP procedure for MS analysis

Sample preparation

- I. It's highly recommended you to provide **IP beads only** and elute protein by our lab. Aspirate extra solution of IP sample.
- II. The recommended starting amount of protein is $\geq 500\mu\text{g}$.
- III. If IP samples are already eluted by yourself, the substance avoid:
SDS, LDS, NP-40, Triton, CHAPS, Octyl glucoside, Octyl thioglucoside, sodium deoxycholate, lauryl maltoside, Brij-35, glycerol, urea, Tris, etc.
- IV. Storage the IP sample lower than -20°C .
- V. You can provide us IP, IP-eluted sample.

Trouble shooting for IP preparation

| High background | |
|--|---|
| Possible cause | Recommend solution |
| Incomplete washing | <ol style="list-style-type: none"> I. General procedure: <ol style="list-style-type: none"> a). Follow your IP procedure to wash step. b). Resuspend the beads with 1X PBS. Transfer the beads to a NEW tube (very important!) and remove all remaining PBS carefully. Store the beads at -80°C. Avoid freeze/thaw cycles. Wash the beads with 1X PBS (Magnetic beads: 10X bed volume; Agarose beads: 20-50X bed volume). c). Repeat step II at least 2-3 times. II. Use more stringent washing buffer for washing: <ol style="list-style-type: none"> a). Using 1% NP-40 but not 0.1 % sodium deoxycholate. Following wash with PBS to remove the detergent. |
| Antibody non-specific binding | Using highly affinity and specificity of purified antibody (e.g., monoclonal antibody) |
| Too many antibodies | Using less antibody for coupling beads |
| Non-specific binding of proteins to antibody | <ol style="list-style-type: none"> I. Reduce the amount of sample II. Pre-clear the lysate by pre-incubating with the beads |
| Too many competing proteins in sample | Spin the lysate for 30 minutes before adding antibody in order to remove debris, insoluble proteins, membrane proteins, other contamination substance, etc. |
| Interfering substances present in sample | <ol style="list-style-type: none"> I. Lysates containing DTT, β-ME, or other reducing agents which can destroy antibody function II. Extreme pH level and high concentration of detergent may also interference antibody and target protein interaction |
| Immunoglobulin contamination | Using crosslinker reagent to decrease antibody dissociate from beads (e.g., BS3 in Dynabeads system) |
| Keratin contamination | <ol style="list-style-type: none"> I. Cleaning the bench and lab equipment before IP preparation II. Using powder-free gloves and new experimental consumables (e.g., filter tips) |

| Weak Signal or No Signal | |
|----------------------------|---|
| Possible cause | Recommend solution |
| Target protein degradation | I. Using fresh lysates sample II. Make sure protease inhibitors are added (e.g., PMSF) |
| Wash too stringent | Reduce the number of washing or decreasing detergent concentration |
| Low titer of antibody | Using titration experiment to determine optimal concentration of antibody |