



TOOLS Bacterial and Fungal DNA and RNA Extraction Kit

**Simple centrifugation procedures for Bacterial and Fungal DNA &
RNA extraction and purification**

Contents

Introduction	2
Kit Contents	2
Workflow.....	3
Protocol	4

Introduction

The TOOLS Bacterial and Fungal DNA and RNA Extraction Kit is a buffer-based system for bacterial and fungal DNA and RNA extraction and purification. The kit involves simple centrifugation procedures that facilitate the complete removal of contaminants and enzyme inhibitors. It affords rapid, simple, and cost-effective extraction and purification processes, and the purified DNA and RNA are suitable for downstream applications, such as polymerase chain reaction, blotting, and sequencing.

Materials not provided

1. 1-Bromo-3-chloropropane (BCP) (CAS Number: 109-70-6).
2. TOOLS Proteinase K (Cat No. RTT-BD03) or other compatible reagents and lysozymes.


Kit Contents

Contents	TX-BFT01 (50 preps)
Buffer BTA	30 ml
Buffer BTB	12 ml
Buffer BTC	300 µl
Binding Gel	50 tubes

Storage

The TOOLS Bacterial and Fungal DNA and RNA Extraction Kit can be stored at room temperature for up to 24 months.

Workflow



Protocol

Sample preparation

1. Place 1 mL of bacterial cells in a 2-mL microcentrifuge tube (RNase free).
2. Centrifuge for 5 min at 14,000 rpm, and remove a sufficient amount of supernatant.
3. Add 500 μ L of BTA buffer to the tube and vortex for 60 s to resuspend the pellet thoroughly.
4. Add 5 μ L of proteinase K (20 mg/mL) and 5 μ L of lysozyme (100 mg/mL) to the tube and vortex for 60 s and then incubate at 55 °C for 1 h. Vortex the tube again for 20 s to mix thoroughly.

DNA and RNA extraction

1. Add 5 μ L of proteinase K and 5 μ L of lysozyme to the tube again.
Vortex for 60 s and then incubate at 55 °C for 30 min.
2. Add 200 μ L of BCP and 200 μ L of BTB buffer to the tube, and mix the sample by inverting the tube three times. Transfer the entire solution to the Binding Gel tube (centrifuge at $13000 \times g$ for 30 s before use).
3. Invert the Binding Gel tube three times (do not vortex) and then centrifuge at $12,000\text{--}16,000 \times g$ for 5 min.
4. Transfer the supernatant to a new microcentrifuge tube, and add 5 μ L of BTC buffer and 550 μ L of isopropanol to the tube.
5. Invert the tube three times and incubate at 37 °C for 10 min.
6. Centrifuge at 14,000 rpm for 5 min. Remove the supernatant.
7. Add 600 μ L of 70% ethanol. Centrifuge at 14,000 rpm for 5 min. Remove the supernatant.
8. Air dry the pellet, and rehydrate the pellet with 20–50 μ L of TE buffer or double-distilled water (ddH₂O) (adjust the buffer volume according to the pellet size).

The product is for research only, not for diagnostic or clinical use.