

## RBC Lysis Buffer

Cat no:	Size
TE-RBC500	500 ml

**Storage:** Room temperature

### Protocol

1. Collect fresh blood in EDTA- $\text{Na}^2$  treated collection tubes (or other anticoagulant mixture).
2. Apply up to 300 $\mu\text{l}$  of blood to a 1.5ml microcentrifuge tube. If blood sample is more than 300 $\mu\text{l}$  (up to 1 ml), apply the sample to a sterile 15ml centrifuge tube.
3. Add 3 times the sample volume of RBC Lysis Buffer and mix by inversion. Do not vortex.
4. Incubate the tube for 5 minutes at room temperature.
5. Centrifuge for 2 minutes at 3,000 x g and discard the supernatant.
6. Add 200 $\mu\text{l}$  RBC Lysis Buffer to resuspend the cell pellet.
7. Repeat step 5 and proceed to further experiments.

The product is for research only, not for diagnostic and clinical use.

