



## TOOLSpeed Blocking Reagent

**Cat. no. TFU-BL500**

**Storage: 4 °C**

**Product Size: 500 mL**

### Introduction

The TOOLSpeed Blocking Reagent is used for effective blocking when performing incubation or hybridization with primary or secondary antibodies.

### Protocol

1. Rinse both sides of PVDF or NC papers with the TOOLSpeed Blocking Reagent for 1–3 min.
2. Dilute primary antibody with the TOOLSpeed Blocking Reagent and perform hybridization with PVDF or NC papers.
3. Wash the membrane with regular wash buffer.
4. Dilute secondary antibody with the TOOLSpeed Blocking Reagent and perform hybridization with PVDF or NC papers.
5. Wash the membrane with regular wash buffer.
6. Proceed to the development procedures with chemiluminescent substrates.

Note: The product is for research purposes only and not for diagnostic and clinical use.