

## dNTP

**Cat. no.:** TG-dNTP、TG-dNTP-2

**Product Size:** 1 mL

| Cat. no.  | Size | Concentration of each nucleotide |
|-----------|------|----------------------------------|
| TG-dNTP   | 1 ml | 2.5 mM                           |
| TG-dNTP-2 | 1 ml | 10 mM                            |

**Storage:** Store at  $-20^{\circ}\text{C}$ , stable for 1 year from the assay date.

**Purity:**  $\geq 99\%$  by HPLC

## Introduction

dNTP is a premixed ready-to-use solution of dATP, dCTP, dGTP, and dTTP (10 mM can be diluted with ddH<sub>2</sub>O in pH 7.5).

## Protocol

This product is suitable for use in all standard PCR techniques, cDNA synthesis, and primer extension.

The product is for research only; not for diagnostic or clinical use.