



## 100 mM dNTP Set (A, T, C, G: 1 mL each)

**Cat. No:** TT-dNTP-100

**Size:** 1 mL each

**Purity:**  $\geq 99.5\%$  purity by HPLC

**Storage:** Stable for 1 year if stored at  $-20^{\circ}\text{C}$ . Avoid multiple freeze/thaw cycles.

## Introduction

dNTP Set contains 0.5 mL of dATP, dCTP, dGTP, and dTTP (100 mM each in sterile deionized water at pH 7.5).

This product is DNase- and RNase-free and suitable for all standard PCR techniques, including DNA sequencing, cDNA synthesis, primer extension, and nick translation.

The product is for research only; not for diagnostic or clinical use.