



TOOLS Tris-Glycine-SDS 10x Buffer

| Cat no: | Size |
|------------|------|
| TEB- BTGS2 | 2 L |

Storage: Room temperature

Description

TOOLS Tris-Glycine-SDS is a biotechnology grade 10X Solution. A single strength (1X) solution contains 0.025 Tris base, 0.192M glycine and 0.1% sodium dodecyl sulfate.

The product is for research only, not for diagnostic and clinical use.

www.tools-biotech.com
info@tools-biotech.com
Your best friend for research

