



#### **4X Upper Buffer**

Cat no:	Size
TEB- BUB01	500 ml

**Storage:** Room temperature

#### **Description**

TOOLS 4X Upper Buffer consist with 0.5M Tris (pH 6.8) and 0.4% SDS at the optimal concentration for convenient and quick preparation of polyacrylamide gels for protein analysis by SDS-PAGE. Dilution of 4X buffer to a 1X solution (0.125M Tris and 0.1% SDS) can be used in the preparation of the stacking gel for SDS-PAGE.

The product is for research only, not for diagnostic and clinical use.

***www.tools-biotech.com***  
***info@tools-biotech.com***  
***Your best friend for research***

