

Test Report

Applicant: Biotools Co., Ltd.
7F-3, No. 99, Sec. 1, Xintai 5th Rd.,
Xizhi District, New Taipei City 221,
Taiwan, R.O.C.

Number : TWNC00660206

Date : Jan 03, 2018

Sample Description:

One (1) group of submitted samples said to be :

Sample Description : DNA View
Style / Item No. : TT-DNA01
Country of Origin : Taiwan
Supplier's Name : Biotools Co., Ltd.
Date Sample Received : Dec 20, 2017
Date Test Started : Dec 20, 2017

Test Conducted :

As requested by the applicant, for details please refer to attached pages.

Authorized by:
On Behalf of Intertek Testing Services
Taiwan Limited



Matt Wang
Sr. Manager



Test Conducted :

1. Allergic disperse dyestuffs(致過敏之分散性染料)

With reference to DIN 54231: 2005, by Liquid Chromatography /Tandem Mass Spectrometer (LC-MS-MS) analysis.

<u>Compound</u>	<u>Result (mg/l In Extract)</u>
	<u>Dark red liquid</u>
Disperse Blue 1	ND
Disperse Blue 3	ND
Disperse Blue 35	ND
Disperse Blue 106	ND
Disperse Blue 124	ND
Disperse Yellow 3	ND
Disperse Orange 3	ND
Disperse Orange 37/76/59	ND
Disperse Red 1	ND

Remarks : Reporting limit = 1 mg/l in extract
ND = Not detected

2. Carcinogenic dyestuffs(致癌性染料)

With reference to DIN 54231: 2005, by Liquid Chromatography /Tandem Mass Spectrometer (LC-MS-MS)analysis.

<u>Compound</u>	<u>Result (mg/l in extract)</u>
	<u>Dark red liquid</u>
Disperse Blue 1	ND
Basic Red 9	ND
Acid Red 26	ND
Disperse Yellow 3	ND
Direct Red 28	ND
Direct Blue 6	ND
Direct Black 38	ND
Disperse Orange 11	ND
Basic Violet 14	ND

Remarks: ND = Not detected
Detection limit = 1 mg/l in extract





End of Report

Except where explicitly agreed in writing, all work and services performed by Intertek is subject to our standard Terms and Conditions which can be obtained at our website: <http://www.intertek-twn.com/terms/>. Intertek's responsibility and liability are limited to the terms and conditions of the agreement.

This report is made solely on the basis of your instructions and / or information and materials supplied by you and provide no warranty on the tested sample(s) be truly representative of the sample source. The report is not intended to be a recommendation for any particular course of action, you are responsible for acting as you see fit on the basis of the report results. Intertek is under no obligation to refer to or report upon any facts or circumstances which are outside the specific instructions received and accepts no responsibility to any parties whatsoever, following the issue of the report, for any matters arising outside the agreed scope of the works. This report does not discharge or release you from your legal obligations and duties to any other person. You are the only one authorized to permit copying or distribution of this report (and then only in its entirety). Any such third parties to whom this report may be circulated rely on the content of the report solely at their own risk.

