



TOOLS SSC 20X buffer

Cat. no.: TEB-20SSC

Storage: Store at room temperature

Product Size: 5L (20X)

Description

TOOLS SSC 20X buffer is saline-sodium citrate (SSC) buffer which is used as a hybridization buffer, to control stringency for washing steps in protocols for Southern blotting, in situ hybridization, DNA Microarray and Northern blotting.

Components

150 mM sodium chloride,

15 mM sodium citrate

Adjust pH to pH 7.0 when dissolve to 1X

The product is for research use only. It is not for diagnostic or clinical use.