



10x PBS Buffer 1L (500 ml x2) STERILE

Cat. no.: TP-PBSS1L

Storage: Store at room temperature

Product Size: 1L

Description

TOOLS 10X PBS is a pH-adjusted phosphate buffers and saline solutions commonly used in biological research, it is diluted to a 1X working concentration (pH should be adjusted).

Components

37mM NaCl, 2.7mM KCl, 8 mM Na₂HPO₄, and 2 mM KH₂PO₄

The product is for research use only. It is not for diagnostic or clinical use.