



### **10x PBS Buffer 1L (500 ml x2) STERILE**

Cat no:	Size
TP-PBSS1L	1L

**Storage:** Room temperature

### **Description**

TOOLS 10X PBS is a pH-adjusted phosphate buffers and saline solutions commonly used in biological research, it is diluted to a 1X working concentration (pH should be adjusted)

### **Components**

37mM NaCl, 2.7mM KCl, 8 mM Na<sub>2</sub>HPO<sub>4</sub>, and 2 mM KH<sub>2</sub>PO<sub>4</sub>

The product is for research only, not for diagnostic and clinical use.

***www.tools-biotech.com***  
***info@tools-biotech.com***  
***Your best friend for research***

