



MSDS (Material Safety Data Sheet)

MSDS revision date: 12/11/2017

1. Product and company identification

Product Name: TOOLS 5X EvaGreen qPCR Mix without ROX

Catalogue Numbers: TF-5EVNR250

Contents

| Product | Content | Cat no: | Size |
|---|---|-------------|------|
| TOOLS 5X EvaGreen qPCR Mix without ROX | TOOLS 5X EvaGreen qPCR Mix with ROX (TOOLS HotTaq DNA polymerase, 5 x qPCR buffer E, 12.5 mM MgCl ₂ (1 x PCR solution – 2.5 mM MgCl ₂), dNTPs, including dTTP to improve reaction sensitivity and efficiency compared to dUTP, EvaGreen® dye) | TF-5EVNR250 | 1 ml |

Class Identifier: For research use only reagent

Company Identification: Biotools Co., Ltd.

7F.-3, No.99, Sec. 1, Xintai 5th Rd., Xizhi Dist., New Taipei City 221, Taiwan (R.O.C.)

Tel: +886-2-26972697 **Fax:** +886-2-26993689 **email:** info@toolsbiotech.com **web:**

<http://www.tools-biotech.com>

2. Composition / information on hazardous ingredients

| Chemical name | % Conc. | CAS# | EC# | Classification (pure ingredient) | Class.(comp.) |
|-------------------------|------------|------|-----|-------------------------------------|---------------|
| Proprietary ingredients | - | - | - | Not hazardous | NA |
| Water | >50% | - | - | Not hazardous | NA |

Chemical characterization

Description: Solution, consisting of the following components.

Other components: Components not listed here are not hazardous or their concentrations do not exceed the limits specified in the OSHA Hazard Communication Standard 29.CFR 1910.1200.

3. Hazards identification

| | |
|------------|--|
| Eyes | May cause eye irritation with susceptible persons. |
| Skin | May cause skin irritation and/or dermatitis. |
| Inhalation | May be harmful by inhalation. |

www.tools-biotech.com
info@tools-biotech.com
Your best friend for research





| | | | | | | | |
|--|---|--------|---|------|---|------------|---|
| Ingestion | May be harmful if swallowed. | | | | | | |
| Carcinogenic effects | No information available. | | | | | | |
| Mutagenic Effects | No information available. | | | | | | |
| Reproductive Toxicity | No information available. | | | | | | |
| Sensitization | No information available. | | | | | | |
| Causes skin and eye irritation. May be harmful by inhalation or swallowed. Health 0 Flammability 0 | | | | | | | |
| Reactivity 0 | | | | | | | |
| HMIS ratings (scale 0 4) | <table><tr><td>HEALTH</td><td>0</td></tr><tr><td>FIRE</td><td>0</td></tr><tr><td>REACTIVITY</td><td>0</td></tr></table> | HEALTH | 0 | FIRE | 0 | REACTIVITY | 0 |
| HEALTH | 0 | | | | | | |
| FIRE | 0 | | | | | | |
| REACTIVITY | 0 | | | | | | |

4. First aid measures

| | |
|--------------------|--|
| Ingestion | Never give anything by mouth to an unconscious person. Consult a physician. |
| After inhalation | Move to fresh air. If symptoms persist, call a physician. |
| After skin contact | Wash off immediately with plenty of water. If symptoms persist, call a physician. |
| After eye contact | Rinse immediately with plenty of water, also under the eyelids, for at least 15 minutes. If eye irritation persists, consult a specialist. |
| Notes to physician | Treat symptomatically. |

5. Fire fighting measures

| | |
|-------------------------------|--|
| Suitable extinguishing agents | Dry chemical. |
| Protective equipment | Wear self-contained breathing apparatus and protective |

6. Accidental release measures

| | |
|-----------------------------------|--|
| Person related safety precautions | Use personal protective equipment. |
| Measures for cleaning/collecting | Soak up with inert absorbent material. |
| Additional information | No dangerous substances are released. |

7. Handling and storage

| | |
|----------|---|
| Handling | Wash thoroughly after handling. Wash hands before eating. Keep container tightly closed. Use with adequate ventilation. |
| Storage | Keep container closed when not in use. Store at room temperature. |

8. Exposure controls and personal protection

| | |
|---|--|
| Components with limit values that require monitoring at the workplace | The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace. |
| Additional information | Ensure adequate ventilation, especially in confined areas. |
| Personal protective equipment | |

www.tools-biotech.com
info@tools-biotech.com
Your best friend for research





| | |
|---------------------------------|--|
| Respiratory Protection | In case of insufficient ventilation wear suitable respiratory equipment. |
| Hand Protection | Protective gloves. |
| Eye Protection | Safety glasses with side-shields. |
| Skin and Body Protection | Lightweight protective clothing. |
| Hygiene Measures | Handle in accordance with good industrial hygiene and safety practice. |
| Environmental Exposure Controls | Prevent product from entering drains. |

9. Physical and chemical protection

| | |
|--------------------------|--|
| General Information | Form: Liquid (buffer); Color: Colorless (buffer); Odor: Odorless |
| Boiling Point/Range | Not applicable. |
| Melting Point/Range | Not applicable. |
| Flash point | Not applicable. |
| Autoignition Temperature | Do not auto ignition. |
| Oxidizing Properties | No information available. |
| Water Solubility | Moderate. |

10. Physical and chemical protection

| | |
|----------------------------------|--|
| Stability | Stable under normal conditions. |
| Materials to be avoided | Metals. |
| Hazardous decomposition products | No information available. |
| Polymerization | Hazardous polymerisation does not occur. |

11. Toxicological information

| | |
|-------------------------|-------------------------------|
| Acute toxicity | No irritant effect. |
| Primary irritant effect | No irritant effect. |
| On the skin | No irritant effect. |
| On the eye | No irritant effect. |
| Sensitization | No sensitizing effects known. |

12. ECOLOGICAL INFORMATION

| | |
|---------------------|---------------------------|
| Ecotoxicity effects | No information available. |
| Mobility | No information available. |
| Biodegradation | No information available. |

www.tools-biotech.com
info@tools-biotech.com
Your best friend for research





Bioaccumulation

No information available.

13. DISPOSAL CONSIDERATIONS

Recommendation Dispose according to current regional and national rules.

14. TRANSPORT INFORMATION

| | |
|------------------------------|---|
| Proper Shipping Name | Not classified as dangerous in the meaning of transport regulations |
| Road/Railway Haulage ADR/RID | Not restricted. |
| Sea Freight IMO (IMDG) | Not restricted. |
| Air Freight IATA (ICAO) | Not restricted. |
| UN Number | Not restricted. |

15. REGULATORY INFORMATION

None of the chemicals in this product are listed as Toxic Pollutants under the CWA.

None of the chemicals in this product are considered highly hazardous by OSHA.

16. Other INFORMATION

This material is sold for research and development purposes only. It is not for any human or animal therapeutic or clinical diagnostic use. It is not intended for food, drug, household, agricultural, or cosmetic use. An individual technically qualified to handle potentially hazardous chemicals must supervise the use of this material.

The above information was acquired by diligent search and/or investigation and the recommendations are based on prudent application of professional judgment. The information shall not be taken as being all inclusive and is to be used only as a guide. All materials and mixtures may be present unknown hazards and should be used with caution. **INFORMATION IN THIS MSDS DOES NOT CONSTITUTE A WARRANTY, EXPRESS OR IMPLIED, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE.**

www.tools-biotech.com
info@tools-biotech.com
Your best friend for research

