



TOOLS siRNA Synthesis Report  
(DEMO REPORT)

服務編號	2100000000001		
Product Name	siRNA		
Label	XXX-siRNA		
Sequence (5'-3')	5' GCUC [REDACTED] UUCGG 3' 5' GAAA [REDACTED] AGCUG 3'		
Modification	N	N	Others: N
Size	2x1 OD		
Reconstitution	To prepare 20μM solution in one tube, please add 125 ul DEPC H <sub>2</sub> O		
Product format	Powder	Purification:	HPLC
Storage	6 months at -20°C		

Important Note:

Synthesized sequence is same as the sequence showed on the synthesis report, if the sequence is different from the one that you provided, please contact with us.

BIOTOOLS Company TEL : +886 2 26963699 e-mail: [info@toolsbiotech.com](mailto:info@toolsbiotech.com)